

CAB20671

Product Information

Product SKU:	CAB20671	Gene ID:	4325329	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Oryza sativa

Additional Information

Observed MW:	100kDa	Conjugate:	Unconjugated
Calculated MW:	127kDa	Isotype:	IgG

Immunogen Information

Background: Involved in regulation of circadian rhythm, and in the control of the photoperiodic flowering. Acts as suppressor of flowering under short-day (SD) and long-day (LD) conditions. Activates Hd1/CONSTANS gene.

Recommended Dilution:	WB, 1:500 - 1:2000
Synonyms:	GI; OsGI
Purification Method:	Affinity purification
Immunogen:	Recombinant protein of Oryza sativa OsGI.
Storage:	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.