

OsWRKY53 Polyclonal Antibody

CAB20672



Product Information

Product SKU:	CAB20672	Gene ID:	4338474	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Oryza sativa

Additional Information

Observed MW:	60kDa	Conjugate:	Unconjugated
Calculated MW:	51kDa	Isotype:	IgG

Immunogen Information

Background: OsWRKY53 is a transcription factor that plays important roles in elicitor-induced defense signaling pathways in rice.

Recommended Dilution:	WB, 1:500 - 1:2000
Synonyms:	OsWRKY53
Purification Method:	Affinity purification
Immunogen:	Recombinant protein of Oryza sativa OsWRKY53.
Storage:	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.