

Product Information

Product SKU:	CAB20675	Gene ID:	4340185	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Oryza sativa

Additional Information

Observed MW:	24kDa	Conjugate:	Unconjugated
Calculated MW:	20kDa	Isotype:	IgG

Immunogen Information

Background: Probable mobile flower-promoting signal (florigen) that moves from the leaf to the shoot apical meristem (SAM) and induces flowering. Promotes the transition from vegetative growth to flowering downstream of HD1 and EHD1 under short day (SD) conditions. Acts upstream of MADS14 and MADS15.

Recommended Dilution:	WB, 1:500 - 1:2000
Synonyms:	Hd3a; OsJ_20191
Purification Method:	Affinity purification
Immunogen:	Recombinant protein of Oryza sativa Hd3a.
Storage:	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.