

Product Information

Product SKU:	CAB20678	Gene ID:	4333000	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Oryza sativa

Additional Information

Observed MW:	26kDa	Conjugate:	Unconjugated
Calculated MW:	26kDa	Isotype:	IgG

Immunogen Information

**Background:** The proteasome is a multicatalytic proteinase complex which is characterized by its ability to cleave peptides with Arg, Phe, Tyr, Leu, and Glu adjacent to the leaving group at neutral or slightly basic pH. The proteasome has an ATP-dependent proteolytic activity.

Recommended Dilution:	WB, 1:500 - 1:2000
Synonyms:	PAB1; OsPAB1; OsJ_11088; OsTT1
Purification Method:	Affinity purification
Immunogen:	Recombinant protein of Oryza sativa OsTT1.
Storage:	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.