

Product Information

Product SKU:	CAB20724	Gene ID:	57466	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	126kDa	Conjugate:	Unconjugated
Calculated MW:	126kDa	Isotype:	IgG

Immunogen Information

Background: This gene likely encodes a member of the arginine/serine-rich splicing factor family. A similar protein in Rat appears to bind the large subunit of RNA polymerase II and provide a link between transcription and pre-mRNA splicing. Alternatively spliced transcript variants have been described.

Recommended Dilution:	IHC-P,1:50 - 1:200
Synonyms:	SRA4; SFRS15; SCAF4
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 170-248 of human SCAF4 (NP_065757.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.