

CAB2092

Product Information

Product SKU: CAB2092	Gene ID: 3934	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human, Mouse

Additional Information

Observed MW: 20kDa	Conjugate: Unconjugated
Calculated MW: 23kDa	Isotype: IgG

Immunogen Information

Background:	This gene encodes a protein that belongs to the lipocalin family. Members of this family transport small hydrophobic molecules such as lipids, steroid hormones and retinoids. The protein encoded by this gene is a neutrophil gelatinase-associated lipocalin and plays a role in innate immunity by limiting bacterial growth as a result of sequestering iron-containing siderophores. The presence of this protein in blood and urine is an early biomarker of acute kidney injury. This protein is thought to be involved in multiple cellular processes, including maintenance of skin homeostasis, and suppression of invasiveness and metastasis. Mice lacking this gene are more susceptible to bacterial infection than wild type mice.
Recommended Dilution:	WB, 1:500 - 1:1000 IHC-P, 1:500 - 1:1000 IF/ICC, 1:50 - 1:200
Synonyms:	p25; 24p3; MSFI; NGAL; Lipocalin-2/NGAL
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 21-198 of human Lipocalin-2/NGAL (NP_005555.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.09% Sodium azide, 50% glycerol, pH7.3.