

Product Information

Product SKU:	CAB20952	Gene ID:	54863	Size:	20uL, 100uL
Clone No:	ARC2939	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	47kDa	Conjugate:	Unconjugated
Calculated MW:	47kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable ATP binding activity. Predicted to be located in extracellular region and platelet alpha granule lumen. Predicted to be active in endoplasmic reticulum lumen and nuclear envelope.

Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	C9orf167; TOR4A
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human TOR4A (Q9NXH8).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.