

CAB21056

Product Information

Product SKU:	CAB21056	Gene ID:	345275	Size:	20uL, 100uL
Clone No:	ARC52038	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	34kDa	Conjugate:	Unconjugated
Calculated MW:	34kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable oxidoreductase activity, acting on the CH-OH group of donors, NAD or NADP as acceptor and steroid dehydrogenase activity. Acts upstream of or within positive regulation of lipid biosynthetic process. Located in lipid droplet.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: FLDP; SCDR9; NIIL497; SDR16C3; HMFN0376; HSD17B13

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 30-211 of human HSD17B13 (NP_835236.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.