
Product Information

Product SKU:	CAB21071	Gene ID:	10023	Size:	20uL, 100uL
Clone No:	ARC2842	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	Refer to figures	Conjugate:	Unconjugated
Calculated MW:	29kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene belongs to the GSK-3-binding protein family. The protein inhibits GSK-3-mediated phosphorylation of beta-catenin and positively regulates the Wnt signaling pathway. It may function in tumor progression and in lymphomagenesis.

Recommended Dilution: IF/ICC, 1:50 - 1:200

Synonyms: FRAT1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 70-170 of human FRAT1 (NP_005470.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.