

---

## Product Information

<b>Product SKU:</b>	CAB21089	<b>Gene ID:</b>	16181	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse

---

## Additional Information

<b>Observed MW:</b>	18kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	20kDa	<b>Isotype:</b>	IgG

---

## Immunogen Information

**Background:** Predicted to enable cytokine receptor binding activity. Acts upstream of or within insulin secretion and lipid metabolic process. Located in extracellular region and vesicle. Human ortholog(s) of this gene implicated in several diseases, including autoimmune disease (multiple); eye disease (multiple); kidney failure (multiple); lung disease (multiple); and vascular disease (multiple). Orthologous to human IL1RN (interleukin 1 receptor antagonist).

**Recommended Dilution:** WB,1:100 - 1:500 IF/ICC,1:50 - 1:200

**Synonyms:** IL-1ra; F630041P17Rik; IL1RN

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 27-178 of mouse IL1RN (NP\_001034790.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.