

CAB21281

Product Information

Product SKU:	CAB21281	Gene ID:	836822	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Arabidopsis thaliana

Additional Information

Observed MW:	41kDa	Conjugate:	Unconjugated
Calculated MW:	41kDa	Isotype:	IgG

Immunogen Information

Background: encodes a member of SNF1-related protein kinases (SnRK2) whose activity is activated by ionic (salt) and non-ionic (mannitol) osmotic stress. Enzyme involved in the ABA signaling during seed germination, dormancy and seedling growth.

Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	SNRK2-3; SRK2I; SUCROSE NONFERMENTING 1 (SNF1)-RELATED PROTEIN KINASE 2-3; sucrose nonfermenting 1(SNF1)-related protein kinase 2.3; SnRK2.3
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 200-300 of arabidopsis thaliana SnRK2.3 (NP_201489.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.