

CAB2151

Product Information

Product SKU:	CAB2151	Gene ID:	64100	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	34kDa	Conjugate:	Unconjugated
Calculated MW:	26kDa	Isotype:	IgG

Immunogen Information

**Background:** The protein encoded by this gene belongs to the sperm-coating protein family of epididymal origin. This protein and its canine homolog are the first known examples of proteins with four tandemly arranged fibronectin type 2 (Fn2) domains in the Fn2-module protein family.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	E12; HE12; EL149; EDDM12; ELSPBP1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-223 of human ELSPBP1 (NP_071425.3).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.