

Product Information

Product SKU:	CAB21681	Gene ID:	84193	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	67kDa	Conjugate:	Unconjugated
Calculated MW:	67kDa	Isotype:	IgG

Immunogen Information

Background: Enables protein-L-histidine N-tele-methyltransferase activity. Involved in actin modification and peptidyl-histidine methylation. Located in cytoplasm.

Recommended Dilution:	WB,1:500 - 1:2000 IP,0.5µg-4µg antibody for 200µg-400µg extracts of whole cells
Synonyms:	hSETD3; C14orf154; SETD3
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 72-329 of human SETD3 (NP_115609.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.