

Cytokeratin 17 (KRT17) Rabbit Polyclonal Antibody

CAB2175



Product Information

Product SKU:	CAB2175	Gene ID:	3872	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	48kDa	Conjugate:	Unconjugated
Calculated MW:	48kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes the type I intermediate filament chain keratin 17, expressed in nail bed, hair follicle, sebaceous glands, and other epidermal appendages. Mutations in this gene lead to Jackson-Lawler type pachyonychia congenita and steatocystoma multiplex.

Recommended Dilution:	WB, 1:500 - 1:2000 IF/ICC, 1:20 - 1:100
Synonyms:	PC; K17; PC2; 39.1; CK-17; PCHC1; Cytokeratin 17 (KRT17)
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 133-432 of human Cytokeratin 17 (Cytokeratin 17 (KRT17)) (NP_000413.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.