

CAB22050

---

## Product Information

<b>Product SKU:</b>	CAB22050	<b>Gene ID:</b>	9975	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC54887	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Mouse,Rat

---

## Additional Information

<b>Observed MW:</b>	70kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	65kDa	<b>Isotype:</b>	IgG

---

## Immunogen Information

**Background:** This gene encodes a member of the nuclear hormone receptor family, specifically the NR1 subfamily of receptors. The encoded protein functions as a transcriptional repressor and may play a role in circadian rhythms and carbohydrate and lipid metabolism. Alternatively spliced transcript variants have been described.

**Recommended Dilution:** WB,1:1000 - 1:5000

**Synonyms:** RVR; BD73; EAR-1R; REVERBB; REVERBbeta; Rev-erb beta/NR1D2

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 220-300 of human Rev-erb beta/NR1D2 (NP\_005117.3).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.