

Product Information

Product SKU:	CAB22117	Gene ID:	282996	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	150kDa	Conjugate:	Unconjugated
Calculated MW:	134kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a protein that binds RNA and regulates splicing. Mutations in this gene have been associated with familial dilated cardiomyopathy.

Recommended Dilution:	WB,1:100 - 1:500
Synonyms:	-
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 820-980 of human RBM20 (NP_001127835.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.