

CAB22262

Product Information

Product SKU:	CAB22262	Gene ID:	79918	Size:	20uL, 100uL
Clone No:	ARC55206	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	53kDa	Conjugate:	Unconjugated
Calculated MW:	53kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a methyltransferase that adds a methyl group to the histone H2AZ, which is involved in nuclear receptor-dependent transcription. The protein also interacts with several endogenous proteins which are involved in nuclear hormone receptor signaling. A related pseudogene is located on chromosome 2. Alternative splicing results in multiple transcript variants.

Recommended Dilution: WB,1:1000 - 1:5000

Synonyms: SETD6

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 310-473 of human SETD6 (NP_001153777.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.