

CAB2341

Product Information

Product SKU:	CAB2341	Gene ID:	-	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	-	Conjugate:	Unconjugated
Calculated MW:	-	Isotype:	IgG

Immunogen Information

Background: 5-Methylcytidine is a modified nucleoside derived from 5-methylcytosine. It is found in ribonucleic acids of animal, plant, and bacterial origin.

Recommended Dilution:	IF/ICC, 1:50 - 1:200
Synonyms:	m5C; N5-Methylcytidine
Purification Method:	Affinity purification
Immunogen:	Chemical compounds corresponding to N5-Methylcytidine / m5C.
Storage:	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.