

TriMethyl-Histone H4-K20 Rabbit Polyclonal Antibody

CAB2372



Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

11kDa

Calculated MW:

11kDa

Applications:

WB IHC IF IP ChIP ChIPseq

Reactivity:

Human, Mouse, Rat, Other
(Wide Range)

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200 IF 1:50 - 1:200 IP
1:50 - 1:200 ChIP 1:20 -
1:100 ChIPseq 1:20 - 1:100

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. This structure consists of approximately 146 bp of DNA wrapped around a nucleosome, an octamer composed of pairs of each of the four core histones (H2A, H2B, H3, and H4). The chromatin fiber is further compacted through the interaction of a linker histone, H1, with the DNA between the nucleosomes to form higher order chromatin structures. This gene is intronless and encodes a replication-dependent histone that is a member of the histone H4 family. Transcripts from this gene lack polyA tails; instead, they contain a palindromic termination element. This gene is found in a histone cluster on chromosome 1. This gene is one of four histone genes in the cluster that are duplicated; this record represents the centromeric copy.

Immunogen information

Gene ID:

8370

Uniprot

P62805

Synonyms:

FO108; H4; H4/n; H4F2; H4FN; HIST2H4; Histone H4; HIST1H4A;
HIST2H4A


Immunogen:

A synthetic methylated peptide corresponding to residues surrounding K20 of human histone H4


Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images




Western blot analysis of extracts of various cell lines, using TriMethyl-Histone H4-K20 antibody (CAB2372) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 1s.



Dot-blot analysis of all sorts of methylation peptides using TriMethyl-Histone H4-K20 antibody (CAB2372).



Immunohistochemistry of paraffin-embedded rat testis using TriMethyl-Histone H4-K20 antibody (CAB2372) at dilution of 1:200 (40x lens).



Immunohistochemistry of paraffin-embedded mouse testis using TriMethyl-Histone H4-K20 antibody (CAB2372) at dilution of 1:200 (40x lens).