

TCTP Rabbit Monoclonal Antibody



CAB2394

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

23kDa

Calculated MW:

20kDa

Applications:

WB IHC IP

Reactivity:

Human, Mouse, Rat

Protein Background

Immunogen information

Gene ID:

7178

Uniprot

P13693

Synonyms:

HRF; TCTP; p02; p23

Immunogen:

A synthesized peptide derived from human TCTP

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC
1:50 - 1:200 IP 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG


Purification:

Affinity purification

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH7.3.

Product Images



Western blot - TCTP Rabbit mAb (CAB2394)