

CAB24518

Product Information

Product SKU:	CAB24518	Gene ID:	9730	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	170kDa	Conjugate:	Unconjugated
Calculated MW:	118kDa/168kDa/169kDa	Isotype:	IgG

Immunogen Information

Background: Enables estrogen receptor binding activity and histone kinase activity (H2A-T120 specific). Involved in cell competition in a multicellular organism; histone H2A-T120 phosphorylation; and negative regulation of transcription by RNA polymerase II. Located in fibrillar center and nucleoplasm. Part of Cul4-RING E3 ubiquitin ligase complex. Colocalizes with COP9 signalosome.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: RIP; VPRBP

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1058-1396 of human VPRBP (NP_055518.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.