

CAB2563

---

## Product Information

<b>Product SKU:</b>	CAB2563	<b>Gene ID:</b>	10971	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

## Additional Information

<b>Observed MW:</b>	32kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	28kDa	<b>Isotype:</b>	IgG

---

## Immunogen Information

**Background:** This gene product belongs to the 14-3-3 family of proteins which mediate signal transduction by binding to phosphoserine-containing proteins. This highly conserved protein family is found in both plants and mammals, and this protein is 99% identical to the mouse and rat orthologs. This gene is upregulated in patients with amyotrophic lateral sclerosis. It contains in its 5' UTR a 6 bp tandem repeat sequence which is polymorphic, however, there is no correlation between the repeat number and the disease.

**Recommended Dilution:** WB,1:500 - 1:2000

**Synonyms:** 1C5; HS1; 14-3-3; 14-3-3 Theta

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide corresponding to a sequence within amino acids 49-149 of human 14-3-3 Theta (NP\_006817.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.