

CAB2611

Product Information

Product SKU:	CAB2611	Gene ID:	2342	Size:	20uL, 100uL
Clone No:	ARC1923	Host Species:	Rabbit	Reactivity:	Human,Rat

Additional Information

Observed MW:	49kDa	Conjugate:	Unconjugated
Calculated MW:	49kDa	Isotype:	IgG

Immunogen Information

**Background:** Enables zinc ion binding activity. Contributes to protein farnesyltransferase activity. Involved in protein farnesylation. Part of microtubule associated complex and protein farnesyltransferase complex.

Recommended Dilution:	WB,1:500 - 1:1000 IHC-P,1:50 - 1:200
Synonyms:	FPTB; FNTB
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human FNTB (P49356).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.