

Caspase-9 Rabbit Polyclonal Antibody



CAB2636

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

36kDa/46kDa

Calculated MW:

17kDa/30kDa/36kDa/46kDa

Applications:

WB IHC IF

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200 IF 1:50 - 1:100

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

This gene encodes a member of the cysteine-aspartic acid protease (caspase) family. Sequential activation of caspases plays a central role in the execution-phase of cell apoptosis. Caspases exist as inactive proenzymes which undergo proteolytic processing at conserved aspartic residues to produce two subunits, large and small, that dimerize to form the active enzyme. This protein can undergo autoproteolytic processing and activation by the apoptosome, a protein complex of cytochrome c and the apoptotic peptidase activating factor 1; this step is thought to be one of the earliest in the caspase activation cascade. This protein is thought to play a central role in apoptosis and to be a tumor suppressor. Alternative splicing results in multiple transcript variants.

Immunogen information

Gene ID:

842

Uniprot

P55211

Synonyms:

CASP9; APAF-3; APAF3; ICE-LAP6; MCH6; PPP1R56; caspase-9;
casp9; Caspase 9

Immunogen:

Recombinant protein of human Caspase-9

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

