

Product Information

Product SKU:	CAB2698	Gene ID:	1260	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	76kDa	Conjugate:	Unconjugated
Calculated MW:	76kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene represents the alpha subunit of a cyclic nucleotide-gated olfactory channel. The encoded protein contains a carboxy-terminal leucine zipper that mediates channel formation.

Recommended Dilution:	WB,1:100 - 1:500
Synonyms:	CNCA; CNG2; CNCA1; OCNC1; OCNCa; OCNALPHA; CNGA2
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CNGA2 (NP_005131.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.