

FOXN2 Rabbit Polyclonal Antibody



CAB2933

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Calculated MW:

37kDa/47kDa

Applications:

IF

Reactivity:

Human

Antibody Information

Recommended dilutions:

IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

This gene encodes a forkhead domain binding protein and may function in the transcriptional regulation of the human T-cell leukemia virus long terminal repeat.

Immunogen information

Gene ID:

3344

Uniprot

P32314

Synonyms:

FOXN2; HTLF


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 172-431 of human FOXN2 (NP_002149.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunofluorescence analysis of A549 cells using FOXN2 antibody (CAB2933).