

# FOXO1 Rabbit Polyclonal Antibody



CAB2934

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

85kDa

### Calculated MW:

69kDa

### Applications:

WB

### Reactivity:

Human, Mouse, Rat

## Protein Background

This gene belongs to the forkhead family of transcription factors which are characterized by a distinct forkhead domain. The specific function of this gene has not yet been determined; however, it may play a role in myogenic growth and differentiation. Translocation of this gene with PAX3 has been associated with alveolar rhabdomyosarcoma.

## Immunogen information

### Gene ID:

2308

### Uniprot

Q12778

### Synonyms:

FKH1; FKHR; FOXO1A; FOXO1

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

### Immunogen:

A synthetic peptide of human FOXO1

### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

