

CAB2996

Product Information

Product SKU:	CAB2996	Gene ID:	3902	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	70kDa	Conjugate:	Unconjugated
Calculated MW:	57kDa	Isotype:	IgG

Immunogen Information

Background: Lymphocyte-activation protein 3 belongs to Ig superfamily and contains 4 extracellular Ig-like domains. The LAG3 gene contains 8 exons. The sequence data, exon/intron organization, and chromosomal localization all indicate a close relationship of LAG3 to CD4.

Recommended Dilution:	WB,1:100 - 1:500
Synonyms:	CD223; LAG3
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 29-170 of human LAG3 (NP_002277.4).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.