

NCAN Rabbit Polyclonal Antibody



CAB3052

Product Information

Size:

50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

143kDa

Applications:

WB IHC

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:25
- 1:100

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

1463

Uniprot

O14594

Synonyms:

NCAN; CSPG3; neurocan

Immunogen:

A synthetic peptide of human NCAN

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
