

CAB3129

Product Information

Product SKU:	CAB3129	Gene ID:	140679	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	53kDa	Conjugate:	Unconjugated
Calculated MW:	57kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene is an integral membrane protein involved in gamma-aminobutyric acid (GABA) and glycine uptake into synaptic vesicles. The encoded protein is a member of amino acid/polyamine transporter family II.

Recommended Dilution: WB,1:1000 - 1:2000

Synonyms: VGAT; VIAAT; SLC32A1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-525 of human SLC32A1 (Q9H598).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.