

PIWIL4 Rabbit Polyclonal Antibody



CAB3144

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

110kDa

Calculated MW:

52kDa/56kDa/96kDa

Applications:

WB IHC

Reactivity:

Human, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

143689

Uniprot

Q7Z3Z4

Synonyms:

PIWIL4; HIWI2; MIWI2

Immunogen:

Recombinant protein of human PIWIL4

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
