

CAB3360

Product Information

Product SKU:	CAB3360	Gene ID:	9692	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Rat

Additional Information

Observed MW:	67kDa	Conjugate:	Unconjugated
Calculated MW:	67kDa	Isotype:	IgG

Immunogen Information

Background: Enables ribonuclease P activity. Involved in mitochondrial tRNA 5'-end processing. Located in mitochondrion and nucleoplasm. Part of mitochondrial ribonuclease P complex.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: MRPP3; COXPD54; KIAA0391

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 314-583 of human KIAA0391 (NP_055487.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.