

# RPA70 Rabbit Monoclonal Antibody



CAB3367

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

70KDa

### Calculated MW:

70kDa

### Applications:

WB

### Reactivity:

Human, Rat

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

6117

#### Uniprot

P27694

#### Synonyms:

HSSB; MST075; REPA1; RF-A; RP-A; RPA70

#### Immunogen:


A synthesized peptide derived from human RPA70

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH7.3.

## Product Images

---



Western blot - RPA70 Rabbit mAb (CAB3367)