

U2SURP Rabbit Polyclonal Antibody

CAB3392



Product Information

Product SKU:	CAB3392	Gene ID:	23350	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	140kDa	Conjugate:	Unconjugated
Calculated MW:	118kDa	Isotype:	IgG

Immunogen Information

Background: Enables RNA binding activity. Predicted to be involved in RNA processing. Located in nucleoplasm.

Recommended Dilution:	WB,1:500 - 1:2000 IHC-P,1:20 - 1:200
Synonyms:	SR140; fSAPa; U2SURP
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 705-1029 of human U2SURP (NP_001073884.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.