

# U2SURP Rabbit Polyclonal Antibody



CAB3392

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

140kDa

### Calculated MW:

72kDa/118kDa

### Applications:

WB IHC IF

### Reactivity:

Human, Mouse, Rat

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:20  
- 1:200 IF 1:50 - 1:200

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

23350

#### Uniprot

O15042

#### Synonyms:

U2SURP; SR140; fSAPa

#### Immunogen:


Recombinant fusion protein containing a sequence corresponding to amino acids 705-1029 of human U2SURP (NP\_001073884.1).

#### Storage:


Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images


---




Western blot analysis of extracts of various cell lines, using U2SURP Antibody (CAB3392) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 15s.



Immunohistochemistry of paraffin-embedded rat liver using U2SURP antibody (CAB3392) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human liver damage using U2SURP antibody (CAB3392) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse liver using U2SURP antibody (CAB3392) at dilution of 1:100 (40x lens).