

Product Information

Product SKU:	CAB3431	Gene ID:	3646	Size:	20uL, 100uL
Clone No:	ARC1997	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	48kDa	Conjugate:	Unconjugated
Calculated MW:	52kDa	Isotype:	IgG

Immunogen Information

**Background:** Enables protein N-terminus binding activity. Contributes to translation initiation factor activity. Involved in positive regulation of mRNA binding activity; regulation of gene expression; and translational initiation. Located in cytosol and nucleus. Part of eukaryotic translation initiation factor 3 complex. Colocalizes with PML body.

**Recommended Dilution:** WB,1:500 - 1:2000 IHC-P,1:50 - 1:200

**Synonyms:** INT6; EIF3S6; EIF3-P48; eIF3-p46; eIF3e

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 210-350 of human eIF3e (P60228).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.