

Product Information

Product SKU:	CAB3432	Gene ID:	51454	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	34kDa	Conjugate:	Unconjugated
Calculated MW:	34kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene is an adapter protein necessary for the engulfment of apoptotic cells by phagocytes. Several transcript variants, some protein coding and some thought not to be protein coding, have been found for this gene.

Recommended Dilution:	WB,1:500 - 1:2000 IHC-P,1:50 - 1:200
Synonyms:	CED6; GULP; CED-6; GULP1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human GULP1 (NP_057399.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.