

Product Information

Product SKU:	CAB3489	Gene ID:	10541	Size:	20uL, 100uL
Clone No:	ARC2014	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	28-31kDa	Conjugate:	Unconjugated
Calculated MW:	29kDa	Isotype:	IgG

Immunogen Information

Background:	Enables RNA polymerase binding activity and histone binding activity. Involved in several processes, including activation of cysteine-type endopeptidase activity involved in apoptotic process; nucleosome assembly; and positive regulation of protein export from nucleus. Located in cytoplasm and nucleoplasm. Colocalizes with nucleolus.
Recommended Dilution:	WB,1:500 - 1:1000 IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:200
Synonyms:	APRIL; SSP29; PHAPI2; ANP32B
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ANP32B (Q92688).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.