

TLE1 Rabbit Monoclonal Antibody



CAB3528

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

90KDa

Calculated MW:

83kDa

Applications:

WB IHC

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC
1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

7088

Uniprot

Q04724

Synonyms:

ESG; ESG1; GRG1


Immunogen:

A synthesized peptide derived from human TLE1

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH7.3.

Product Images



Western blot - TLE1 Rabbit mAb (CAB3528)