

Product Information

Product SKU:	CAB3543	Gene ID:	126626	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	55kDa	Conjugate:	Unconjugated
Calculated MW:	49kDa	Isotype:	IgG

Immunogen Information

Background: Enables transcription cis-regulatory region binding activity. Acts upstream of or within positive regulation of transcription by RNA polymerase II. Located in nucleus.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	GABPB-2; GABPB2
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 239-448 of human GABPB2 (NP_653219.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.