

CAB3551

Product Information

Product SKU:	CAB3551	Gene ID:	150864	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	65kDa	Conjugate:	Unconjugated
Calculated MW:	62kDa	Isotype:	IgG

Immunogen Information

Background: FAM117B (Family With Sequence Similarity 117 Member B) is a Protein Coding gene. Diseases associated with FAM117B include Lateral Sclerosis and Amyotrophic Lateral Sclerosis 2, Juvenile. An important paralog of this gene is GLCC11.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	ALS2CR13; FAM117B
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 293-372 of human FAM117B (NP_775782.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.