

CAB3597

Product Information

Product SKU:	CAB3597	Gene ID:	668	Size:	20uL, 100uL
Clone No:	ARC2055	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	50kDa	Conjugate:	Unconjugated
Calculated MW:	39kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a forkhead transcription factor. The protein contains a fork-head DNA-binding domain and may play a role in ovarian development and function. Expansion of a polyalanine repeat region and other mutations in this gene are a cause of blepharophimosis syndrome and premature ovarian failure 3.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: BPES; PFRK; POF3; BPES1; PINTO; FOXL2

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human FOXL2 (P58012).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.