

---

## Product Information

<b>Product SKU:</b>	CAB3667	<b>Gene ID:</b>	11031	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

## Additional Information

<b>Observed MW:</b>	22kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	22kDa	<b>Isotype:</b>	IgG

---

## Immunogen Information

**Background:** Enables GDP binding activity and GTP binding activity. Involved in several processes, including Golgi to plasma membrane protein transport; cellular response to insulin stimulus; and receptor internalization. Located in early endosome; phagocytic vesicle; and trans-Golgi network membrane. Biomarker of severe acute respiratory syndrome.

**Recommended Dilution:** WB,1:200 - 1:2000 IF/ICC,1:50 - 1:200

**Synonyms:** Rab22B; RAB31

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 1-195 of human RAB31 (NP\_006859.2).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.