

Product Information

Product SKU:	CAB3682	Gene ID:	56605	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	48kDa	Conjugate:	Unconjugated
Calculated MW:	54kDa	Isotype:	IgG

Immunogen Information

Background: Enables thiol oxidase activity. Involved in protein folding in endoplasmic reticulum. Predicted to be located in membrane. Predicted to be active in endoplasmic reticulum membrane.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	ERO1LB; Ero1beta
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 34-230 of human ERO1LB (NP_063944.3).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.