

## CAB3800

---

**Product Information**

<b>Product SKU:</b>	CAB3800	<b>Gene ID:</b>	4841	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC0836	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

**Additional Information**

<b>Observed MW:</b>	60kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	54kDa	<b>Isotype:</b>	IgG

---

**Immunogen Information**

**Background:** This gene encodes an RNA-binding protein which plays various roles in the nucleus, including transcriptional regulation and RNA splicing. A rearrangement between this gene and the transcription factor E3 gene has been observed in papillary renal cell carcinoma. Alternatively spliced transcript variants have been described. Pseudogenes exist on Chromosomes 2 and 16.

**Recommended Dilution:** WB,1:500 - 1:1000 IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:200

**Synonyms:** P54; NMT55; NRB54; MRXS34; P54NRB; PPP1R114; NONO/p54nrb

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide corresponding to a sequence within amino acids 372-471 of human NONO/p54nrb (Q15233).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.