

CAB3808

Product Information

Product SKU:	CAB3808	Gene ID:	10399	Size:	20uL, 100uL
Clone No:	ARC0837	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	32kDa	Conjugate:	Unconjugated
Calculated MW:	35kDa	Isotype:	IgG

Immunogen Information

Background: Enables several functions, including cyclin binding activity; enzyme binding activity; and protein domain specific binding activity. Involved in several processes, including positive regulation of hydrolase activity; regulation of cellular protein metabolic process; and regulation of signal transduction. Located in several cellular components, including midbody; perinuclear region of cytoplasm; and phagocytic cup. Part of IRE1-RACK1-PP2A complex.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: H12.3; HLC-7; PIG21; GNB2L1; Gnb2-rs1; RACK1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 220-310 of human RACK1 (P63244).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.