

RACK1 Rabbit Monoclonal Antibody



CAB3808

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

32KDa

Calculated MW:

35kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

10399

Uniprot

P63244

Synonyms:

GNB2L1; Gnb2-rs1; H12.3; HLC-7; PIG21


Immunogen:

A synthesized peptide derived from human RACK1

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH7.3.

Product Images



Western blot - RACK1 Rabbit mAb (CAB3808)