

Product Information

Product SKU:	CAB3842	Gene ID:	1857	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	93kDa	Conjugate:	Unconjugated
Calculated MW:	78kDa	Isotype:	IgG

Immunogen Information

Background: This gene is a member of a multi-gene family which shares strong similarity with the Drosophila dishevelled gene, dsh. The Drosophila dishevelled gene encodes a cytoplasmic phosphoprotein that regulates cell proliferation.

Recommended Dilution:	WB,1:500 - 1:1000 IHC-P,1:50 - 1:200
Synonyms:	DRS3; DVL3
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 580-650 of human DVL3 (NP_004414.3).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.