

PIK3C3 Rabbit Polyclonal Antibody



CAB4021

Product Information

Size:

100uL, 200uL

Observed MW:

100kDa

Calculated MW:

101kDa

Applications:

WB IHC IF

Reactivity:

Human, Mouse

Antibody Information

Recommended dilutions:

WB 1:200 - 1:2000 IHC 1:50
- 1:200 IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

5289

Uniprot

Q8NEB9

Synonyms:

PIK3C3; VPS34; Vps34; hVps34


Immunogen:

A synthetic peptide of human PIK3C3

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of extracts of various cell lines, using PIK3C3 antibody (CAB4021). Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.