

CAB4091

---

## Product Information

<b>Product SKU:</b>	CAB4091	<b>Gene ID:</b>	6426	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC51453	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

## Additional Information

<b>Observed MW:</b>	33kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	28kDa	<b>Isotype:</b>	IgG

---

## Immunogen Information

**Background:** This gene encodes a member of the arginine/serine-rich splicing factor protein family. The encoded protein can either activate or repress splicing, depending on its phosphorylation state and its interaction partners. Multiple transcript variants have been found for this gene. There is a pseudogene of this gene on chromosome 13.

**Recommended Dilution:** WB,1:1000 - 1:5000 IHC-P,1:500 - 1:1000

**Synonyms:** ASF; SF2; SFRS1; SF2p33; SRp30a; SRSF1/SF2/ASF

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide corresponding to a sequence within amino acids 100-200 of human SRSF1/SF2/ASF (NP\_008855.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.